



RU Healthy?

Technical Documentation

Class

ECE 567 – Software Engineering 1

Semester

Fall 2017

Group #2

Tahiya Chowdhury, Tina Drew,
George Koubbe, Aymen Al-Saadi,
Himabindu Paruchuri and Ramya Tadepalli

Website

[RU Healthy?](#)

Generated by Aymen and
George 1.8.13

Contents

1	Namespace Index	1
1.1	Packages	1
2	Hierarchical Index	3
2.1	Class Hierarchy	3
3	Class Index	5
3.1	Class List	5
4	File Index	7
4.1	File List	7
5	Namespace Documentation	9
5.1	Package com	9
5.2	Package com.example	9
5.3	Package com.example.RUHEALTHY	9
5.4	Package com.example.RUHEALTHY.demo1	9
6	Class Documentation	11
6.1	com.example.RUHEALTHY.demo1.DashboardActivity Class Reference	11
6.1.1	Member Function Documentation	11
6.1.1.1	onCreate()	11
6.1.2	Member Data Documentation	11
6.1.2.1	Email	12
6.1.2.2	EmailHolder	12
6.1.2.3	HRM	12

6.1.2.4	LogOUT	12
6.1.2.5	Profile	12
6.1.2.6	Step	12
6.2	com.example.RUHEALTHY.demo1.HeartbeatView Class Reference	12
6.2.1	Detailed Description	13
6.2.2	Constructor & Destructor Documentation.....	13
6.2.2.1	HeartbeatView() [1 / 2]	13
6.2.2.2	HeartbeatView() [2 / 2]	13
6.2.3	Member Function Documentation	13
6.2.3.1	onDraw()	14
6.2.3.2	onMeasure()	14
6.2.4	Member Data Documentation	14
6.2.4.1	greenBitmap	14
6.2.4.2	matrix.....	14
6.2.4.3	paint.....	14
6.2.4.4	parentHeight	14
6.2.4.5	parentWidth	14
6.2.4.6	redBitmap	15
6.3	com.example.RUHEALTHY.demo1.HeartRateMonitor Class Reference	15
6.3.1	Detailed Description	16
6.3.2	Member Function Documentation	16
6.3.2.1	getCurrent()	16
6.3.2.2	getSmallestPreviewSize()	16
6.3.2.3	onConfigurationChanged()	16
6.3.2.4	onCreate()	16
6.3.2.5	onPause()	16
6.3.2.6	onResume()	17
6.3.3	Member Data Documentation	17
6.3.3.1	averageArray	17
6.3.3.2	averageArraySize	17

6.3.3.3	averageIndex	17
6.3.3.4	beats	17
6.3.3.5	beatsArray	17
6.3.3.6	beatsArraySize	17
6.3.3.7	beatsIndex	18
6.3.3.8	camera.....	18
6.3.3.9	currentType.....	18
6.3.3.10	image.....	18
6.3.3.11	preview	18
6.3.3.12	previewCallback.....	18
6.3.3.13	previewHolder.....	18
6.3.3.14	processing	19
6.3.3.15	startTime.....	19
6.3.3.16	surfaceCallback	19
6.3.3.17	TAG	19
6.3.3.18	text.....	19
6.3.3.19	wakeLock.....	19
6.4	com.example.RUHEALTHY.demo1.ImageProcessing Class Reference	19
6.4.1	Detailed Description	20
6.4.2	Member Function Documentation	20
6.4.2.1	decodeYUV420SPtoRedAvg().....	20
6.4.2.2	decodeYUV420SPtoRedSum().....	20
6.5	com.example.RUHEALTHY.demo1.MainActivity Class Reference	21
6.5.1	Member Function Documentation	21
6.5.1.1	CheckEditTextStatus().....	21
6.5.1.2	CheckFinalResult()	22
6.5.1.3	LoginFunction().....	22
6.5.1.4	onCreate()	22
6.5.2	Member Data Documentation	22
6.5.2.1	cursor.....	22

6.5.2.2	EditTextEmptyHolder	22
6.5.2.3	Email	22
6.5.2.4	EmailHolder	22
6.5.2.5	LogInButton	23
6.5.2.6	Password	23
6.5.2.7	PasswordHolder	23
6.5.2.8	RegisterButton	23
6.5.2.9	sqliteDatabaseObj	23
6.5.2.10	sqliteHelper	23
6.5.2.11	TempPassword	23
6.5.2.12	UserEmail	24
6.6	com.example.RUHEALTHY.demo1.ProfileActivity Class Reference	24
6.6.1	Member Function Documentation	24
6.6.1.1	CheckEditTextStatus()	25
6.6.1.2	InsertDataIntoSQLiteDatabase()	25
6.6.1.3	onCreate()	25
6.6.1.4	SQLiteDataBaseBuild()	25
6.6.2	Member Data Documentation	25
6.6.2.1	Age	25
6.6.2.2	AgeHolder	25
6.6.2.3	EditTextEmptyHolder	25
6.6.2.4	EmailHolder	26
6.6.2.5	Gender	26
6.6.2.6	GenderHolder	26
6.6.2.7	Height	26
6.6.2.8	HeightHolder	26
6.6.2.9	sqliteDatabaseObj	26
6.6.2.10	SQLiteDataBaseQueryHolder	26
6.6.2.11	sqliteHelper	26
6.6.2.12	Submit	27

6.6.2.13	Target	27
6.6.2.14	TargetHolder	27
6.6.2.15	Weight	27
6.6.2.16	WeightHolder	27
6.7	com.example.RUHEALTHY.demo1.RegisterActivity Class Reference	27
6.7.1	Member Function Documentation	28
6.7.1.1	CheckEditTextStatus()	28
6.7.1.2	CheckFinalResult()	28
6.7.1.3	CheckingEmailAlreadyExistsOrNot()	28
6.7.1.4	EmptyEditTextAfterDataInsert()	28
6.7.1.5	InsertDataIntoSQLiteDatabase()	29
6.7.1.6	onCreate()	29
6.7.1.7	SQLiteDataBaseBuild()	29
6.7.1.8	SQLiteTableBuild()	29
6.7.2	Member Data Documentation	29
6.7.2.1	cursor	29
6.7.2.2	EditTextEmptyHolder	29
6.7.2.3	Email	29
6.7.2.4	EmailHolder	30
6.7.2.5	F_Result	30
6.7.2.6	Name	30
6.7.2.7	NameHolder	30
6.7.2.8	Password	30
6.7.2.9	PasswordHolder	30
6.7.2.10	Register	30
6.7.2.11	sqliteDatabaseObj	30
6.7.2.12	SQLiteDataBaseQueryHolder	31
6.7.2.13	sqliteHelper	31
6.8	com.example.RUHEALTHY.demo1.SensorFilter Class Reference	31
6.8.1	Detailed Description	31

6.8.2	Constructor & Destructor Documentation.....	31
6.8.2.1	SensorFilter().....	31
6.8.3	Member Function Documentation	31
6.8.3.1	cross().....	32
6.8.3.2	dot()	32
6.8.3.3	norm()	32
6.8.3.4	normalize().....	32
6.8.3.5	sum()	32
6.9	com.example.RUHEALTHY.demo1.SQLiteHelper Class Reference	32
6.9.1	Detailed Description	33
6.9.2	Constructor & Destructor Documentation.....	33
6.9.2.1	SQLiteHelper().....	33
6.9.3	Member Function Documentation	33
6.9.3.1	onCreate()	33
6.9.3.2	onUpgrade().....	34
6.9.4	Member Data Documentation	34
6.9.4.1	DATABASE_NAME	34
6.9.4.2	Table_Column_1_Name.....	34
6.9.4.3	Table_Column_2_Email	34
6.9.4.4	Table_Column_3_Password.....	34
6.9.4.5	Table_Column_4_Age.....	34
6.9.4.6	Table_Column_5_Gender	34
6.9.4.7	Table_Column_6_Height	35
6.9.4.8	Table_Column_7_Weight	35
6.9.4.9	Table_Column_8_Target	35
6.9.4.10	TABLE_NAME	35
6.10	com.example.RUHEALTHY.demo1.StepCounter Class Reference	35
6.10.1	Member Function Documentation	36
6.10.1.1	onAccuracyChanged()	36
6.10.1.2	onCreate()	36

6.10.1.3	onSensorChanged()	36
6.10.1.4	step()	36
6.10.2	Member Data Documentation	37
6.10.2.1	accel	37
6.10.2.2	calories	37
6.10.2.3	miles	37
6.10.2.4	numSteps	37
6.10.2.5	sensorManager	37
6.10.2.6	simpleStepDetector	37
6.10.2.7	TEXT_Miles_Calorie	37
6.10.2.8	TEXT_Miles_STEPS	38
6.10.2.9	TEXT_NUM_STEPS	38
6.10.2.10	TvCalories	38
6.10.2.11	TvMiles	38
6.10.2.12	TvSteps	38
6.11	com.example.RUHEALTHY.demo1.StepDetector Class Reference	38
6.11.1	Detailed Description	39
6.11.2	Member Function Documentation	39
6.11.2.1	registerListener()	39
6.11.2.2	updateAccel()	39
6.11.3	Member Data Documentation	39
6.11.3.1	ACCEL_RING_SIZE	40
6.11.3.2	accelRingCounter	40
6.11.3.3	accelRingX	40
6.11.3.4	accelRingY	40
6.11.3.5	accelRingZ	40
6.11.3.6	lastStepTimeNs	40
6.11.3.7	listener	40
6.11.3.8	oldVelocityEstimate	41
6.11.3.9	STEP_DELAY_NS	41
6.11.3.10	STEP_THRESHOLD	41
6.11.3.11	VEL_RING_SIZE	41
6.11.3.12	velRing	41
6.11.3.13	velRingCounter	41
6.12	com.example.RUHEALTHY.demo1.StepListener Interface Reference	41
6.12.1	Detailed Description	42
6.12.2	Member Function Documentation	42
6.12.2.1	step()	42
6.13	com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE Enum Reference	42
6.13.1	Member Data Documentation	42
6.13.1.1	GREEN	42
6.13.1.2	RED	42

7 File Documentation	43
7.1 DashboardActivity.java File Reference	43
7.2 HeartbeatView.java File Reference	43
7.3 HeartRateMonitor.java File Reference.....	43
7.4 ImageProcessing.java File Reference	44
7.5 MainActivity.java File Reference.....	44
7.6 ProfileActivity.java File Reference	44
7.7 RegisterActivity.java File Reference	44
7.8 SensorFilter.java File Reference.....	45
7.9 SQLiteHelper.java File Reference	45
7.10 StepCounter.java File Reference.....	45
7.11 StepDetector.java File Reference	45
7.12 StepListener.java File Reference.....	45
Index	47

Chapter 1

Namespace Index

1.1 Packages

Here are the packages with brief descriptions (if available):

com	9
com.example	9
com.example.RUHEALTHY	9
com.example.RUHEALTHY.demo1	9

Chapter 2

Hierarchical Index

2.1 Class Hierarchy

com.example.RUHEALTHY.demo1.ImageProcessing	19
com.example.RUHEALTHY.demo1.SensorFilter	31
com.example.RUHEALTHY.demo1.StepDetector	38
com.example.RUHEALTHY.demo1.StepListener	41
com.example.RUHEALTHY.demo1.StepCounter	35
com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE	42
Activity	
com.example.RUHEALTHY.demo1.HeartRateMonitor.	15
AppCompatActivity	
com.example.RUHEALTHY.demo1.DashboardActivity	11
com.example.RUHEALTHY.demo1.MainActivity	21
com.example.RUHEALTHY.demo1.ProfileActivity	24
com.example.RUHEALTHY.demo1.RegisterActivity	27
com.example.RUHEALTHY.demo1.StepCounter	35
SensorEventListener	
com.example.RUHEALTHY.demo1.StepCounter	35
SQLiteOpenHelper	
com.example.RUHEALTHY.demo1.SQLiteHelper	32
..	
View	
com.example.RUHEALTHY.demo1.HeartbeatView	12

Chapter 3

Class Index

3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

com.example.RUHEALTHY.demo1.DashboardActivity	11
com.example.RUHEALTHY.demo1.HeartbeatView	12
com.example.RUHEALTHY.demo1.HeartRateMonitor	15
com.example.RUHEALTHY.demo1.ImageProcessing	19
com.example.RUHEALTHY.demo1.MainActivity	21
com.example.RUHEALTHY.demo1.ProfileActivity	24
com.example.RUHEALTHY.demo1.RegisterActivity	27
com.example.RUHEALTHY.demo1.SensorFilter	31
com.example.RUHEALTHY.demo1.SQLiteHelper	32
com.example.RUHEALTHY.demo1.StepCounter	35
com.example.RUHEALTHY.demo1.StepDetector	38
com.example.RUHEALTHY.demo1.StepListener	41
com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE	42

Chapter 4

File Index

4.1 File List

Here is a list of all files with brief descriptions:

DashboardActivity.java	43
HeartbeatView.java	43
HeartRateMonitor.java	43
ImageProcessing.java	44
MainActivity.java	44
ProfileActivity.java	44
RegisterActivity.java	44
SensorFilter.java	45
SQLiteHelper.java	45
StepCounter.java	45
StepDetector.java	45
StepListener.java	45

Chapter 5

Namespace Documentation

5.1 Package com

Packages

- package **example**

5.2 Packagecom.example

Packages

- package **RUHEALTHY**

5.3 Packagecom.example.RUHEALTHY

Packages

- package **demo1**

5.4 Packagecom.example.RUHEALTHY.demo1

Classes

- class **DashboardActivity**
- class **HeartbeatView**
- class **HeartRateMonitor**
- class **ImageProcessing**
- class **MainActivity**
- class **ProfileActivity**
- class **RegisterActivity**
- class **SensorFilter**
- class **SQLiteHelper**
- class **StepCounter**
- class **StepDetector**
- interface **StepListener**

Chapter 6

Class Documentation

6.1 com.example.RUHEALTHY.demo1.DashboardActivityClass Reference

Inherits AppCompatActivity.

Protected Member Functions

- void **onCreate** (Bundle savedInstanceState)

Package Attributes

- String **EmailHolder**
- TextView **Email**
- Button **LogOUT**
- Button **Profile**
- Button **HRM**
- Button **Step**

6.1.1 Member Function Documentation

6.1.1.1 onCreate()

```
void      com.example.RUHEALTHY.demo1.DashboardActivity.onCreate      (  
          Bundle savedInstanceState ) [protected]
```

6.1.2 Member Data Documentation

6.1.2.1 Email

`TextView com.example.RUHEALTHY.demo1.DashboardActivity.Email [package]`

6.1.2.2 EmailHolder

`String com.example.RUHEALTHY.demo1.DashboardActivity.EmailHolder [package]`

6.1.2.3 HRM

`Button com.example.RUHEALTHY.demo1.DashboardActivity.HRM [package]`

6.1.2.4 LogOUT

`Button com.example.RUHEALTHY.demo1.DashboardActivity.LogOUT [package]`

6.1.2.5 Profile

`Button com.example.RUHEALTHY.demo1.DashboardActivity.Profile [package]`

6.1.2.6 Step

`Button com.example.RUHEALTHY.demo1.DashboardActivity.Step [package]`

The documentation for this class was generated from the following file:

- **DashboardActivity.java**

6.2 com.example.RUHEALTHY.demo1.HeartbeatView ClassReference

Inherits View.

Public Member Functions

- **HeartbeatView** (Context context, AttributeSet attr)
- **HeartbeatView** (Context context)

Protected Member Functions

- void **onMeasure** (int widthMeasureSpec, int heightMeasureSpec)
- void **onDraw** (Canvas canvas)

Static Private Attributes

- static final Matrix **matrix** = new Matrix()
- static final Paint **paint** = new Paint(Paint.ANTI_ALIAS_FLAG)
- static Bitmap **greenBitmap** = null
- static Bitmap **redBitmap** = null
- static int **parentWidth** = 0
- static int **parentHeight** = 0

6.2.1 Detailed Description

This class extends the View class and is designed draw the heartbeat image.

Author

Justin Wetherell phishman3579@gmail.com

6.2.2 Constructor & Destructor Documentation

6.2.2.1 HeartbeatView() [1 / 2]

```
com.example.RUHEALTHY.demo1.HeartbeatView.HeartbeatView (
    Context context,
    AttributeSet attr )
```

6.2.2.2 HeartbeatView() [2 / 2]

```
com.example.RUHEALTHY.demo1.HeartbeatView.HeartbeatView (
    Context context )
```

6.2.3 Member Function Documentation

6.2.3.1 onDraw()

```
void    com.example.RUHEALTHY.demol.HeartbeatView.onDraw    (  
        Canvas canvas ) [protected]
```

6.2.3.2 onMeasure()

```
void    com.example.RUHEALTHY.demol.HeartbeatView.onMeasure  (  
        int widthMeasureSpec,  
        int heightMeasureSpec ) [protected]
```

6.2.4 Member Data Documentation

6.2.4.1 greenBitmap

```
Bitmap com.example.RUHEALTHY.demol.HeartbeatView.greenBitmap = null [static],[private]
```

6.2.4.2 matrix

```
final Matrix com.example.RUHEALTHY.demol.HeartbeatView.matrix = new Matrix() [static],[private]
```

6.2.4.3 paint

```
final Paint com.example.RUHEALTHY.demol.HeartbeatView.paint = new Paint(Paint.ANTI_ALIAS_FLAG)  
[static], [private]
```

6.2.4.4 parentHeight

```
int com.example.RUHEALTHY.demol.HeartbeatView.parentHeight = 0 [static],[private]
```

6.2.4.5 parentWidth

```
int com.example.RUHEALTHY.demol.HeartbeatView.parentWidth = 0 [static],[private]
```

6.2.4.6 redBitmap

```
Bitmap com.example.RUHEALTHY.demo1.HeartbeatView.redBitmap = null [static], [private]
```

The documentation for this class was generated from the following file:

- **HeartbeatView.java**

6.3 com.example.RUHEALTHY.demo1.HeartRateMonitor ClassReference

Inherits Activity.

Collaboration diagram for com.example.RUHEALTHY.demo1.HeartRateMonitor:

Classes

- enum **TYPE**

Public Member Functions

- void **onCreate** (Bundle savedInstanceState)
- void **onConfigurationChanged** (Configuration newConfig)
- void **onResume** ()
- void **onPause** ()

Static Public Member Functions

- static **TYPE** **getCurrent** ()

Static Private Member Functions

- static Camera.Size **getSmallestPreviewSize** (int width, int height, Camera.Parameters parameters)

Static Private Attributes

- static final String **TAG** = "HeartRateMonitor"
- static final AtomicBoolean **processing** = new AtomicBoolean(false)
- static SurfaceView **preview** =null
- static SurfaceHolder **previewHolder** = null
- static Camera **camera** =null
- static View **image** =null
- static TextView **text** = null
- static WakeLock **wakeLock** =null
- static int **averageIndex** =0
- static final int **averageArraySize** =4
- static final int [] **averageArray** = new int[**averageArraySize**]
- static **TYPE** **currentType** = **TYPE.GREEN**
- static int **beatsIndex** =0
- static final int **beatsArraySize** =3
- static final int [] **beatsArray** = new int[**beatsArraySize**]
- static double **beats** =0
- static long **startTime** =0
- static PreviewCallback **previewCallback**
- static SurfaceHolder.Callback **surfaceCallback**

6.3.1 Detailed Description

This class extends Activity to handle a picture preview, process the preview for a red values and determine a heart beat.

Author

Justin Wetherell phishman3579@gmail.com

6.3.2 Member Function Documentation

6.3.2.1 getCurrent()

```
static TYPE com.example.RUHEALTHY.demol.HeartRateMonitor.getCurrent ( ) [static]
```

Here is the caller graph for this function:

6.3.2.2 getSmallestPreviewSize()

```
static Camera.Size com.example.RUHEALTHY.demol.HeartRateMonitor.getSmallestPreviewSize (
    int width,
    int height,
    Camera.Parameters parameters ) [static], [private]
```

6.3.2.3 onConfigurationChanged()

```
void com.example.RUHEALTHY.demol.HeartRateMonitor.onConfigurationChanged (
    Configuration newConfig )
```

6.3.2.4 onCreate()

```
void com.example.RUHEALTHY.demol.HeartRateMonitor.onCreate (
    Bundle savedInstanceState )
```

6.3.2.5 onPause()

```
void com.example.RUHEALTHY.demol.HeartRateMonitor.onPause ( )
```

6.3.2.6 onResume()

```
void com.example.RUHEALTHY.demo1.HeartRateMonitor.onResume ( )
```

6.3.3 Member Data Documentation

6.3.3.1 averageArray

```
final int [ ] com.example.RUHEALTHY.demo1.HeartRateMonitor.averageArray = new int[ average←  
ArraySize] [static], [private]
```

6.3.3.2 averageArraySize

```
final int com.example.RUHEALTHY.demo1.HeartRateMonitor.averageArraySize = 4 [static], [private]
```

6.3.3.3 averageIndex

```
int com.example.RUHEALTHY.demo1.HeartRateMonitor.averageIndex = 0 [static], [private]
```

6.3.3.4 beats

```
double com.example.RUHEALTHY.demo1.HeartRateMonitor.beats = 0 [static], [private]
```

6.3.3.5 beatsArray

```
final int [ ] com.example.RUHEALTHY.demo1.HeartRateMonitor.beatsArray = new int[ beatsArray←  
Size] [static], [private]
```

6.3.3.6 beatsArraySize

```
final int com.example.RUHEALTHY.demo1.HeartRateMonitor.beatsArraySize = 3 [static], [private]
```

6.3.3.7 beatsIndex

```
int com.example.RUHEALTHY.demo1.HeartRateMonitor.beatsIndex = 0 [static],[private]
```

6.3.3.8 camera

```
Camera com.example.RUHEALTHY.demo1.HeartRateMonitor.camera = null [static],[private]
```

6.3.3.9 currentType

```
TYPE com.example.RUHEALTHY.demo1.HeartRateMonitor.currentType = TYPE.GREEN [static],[private]
```

6.3.3.10 image

```
View com.example.RUHEALTHY.demo1.HeartRateMonitor.image = null [static],[private]
```

6.3.3.11 preview

```
SurfaceView com.example.RUHEALTHY.demo1.HeartRateMonitor.preview = null [static],[private]
```

6.3.3.12 previewCallback

```
PreviewCallback com.example.RUHEALTHY.demo1.HeartRateMonitor.previewCallback [static],[private]
```

6.3.3.13 previewHolder

```
SurfaceHolder com.example.RUHEALTHY.demo1.HeartRateMonitor.previewHolder = null [static],[private]
```

6.3.3.14 processing

```
final AtomicBoolean com.example.RUHEALTHY.demo1.HeartRateMonitor.processing = new AtomicBoolean(false) [static], [private]
```

6.3.3.15 startTime

```
long com.example.RUHEALTHY.demo1.HeartRateMonitor.startTime = 0 [static], [private]
```

6.3.3.16 surfaceCallback

```
SurfaceHolder.Callback com.example.RUHEALTHY.demo1.HeartRateMonitor.surfaceCallback [static], [private]
```

6.3.3.17 TAG

```
final String com.example.RUHEALTHY.demo1.HeartRateMonitor.TAG = "HeartRateMonitor" [static], [private]
```

6.3.3.18 text

```
TextView com.example.RUHEALTHY.demo1.HeartRateMonitor.text = null [static], [private]
```

6.3.3.19 wakeLock

```
WakeLock com.example.RUHEALTHY.demo1.HeartRateMonitor.wakeLock = null [static], [private]
```

The documentation for this class was generated from the following file:

- **HeartRateMonitor.java**

6.4 com.example.RUHEALTHY.demo1.ImageProcessing Class Reference

Static Public Member Functions

- static int **decodeYUV420SPtoRedAvg** (byte[] yuv420sp, int width, int height)

Static Private Member Functions

- static int **decodeYUV420SPtoRedSum** (byte[] yuv420sp, int width, int height)

6.4.1 Detailed Description

This abstract class is used to process images.

Author

Justin Wetherell phishman3579@gmail.com

6.4.2 Member Function Documentation

6.4.2.1 decodeYUV420SPtoRedAvg()

```
static int com.example.RUHEALTHY.demol.ImageProcessing.decodeYUV420SPtoRedAvg (
    byte [] yuv420sp,
    int width,
    int height ) [static]
```

Given a byte array representing a yuv420sp image, determine the average amount of red in the image. Note: returns 0 if the byte array is NULL.

Parameters

<i>yuv420sp</i>	Byte array representing a yuv420sp image
<i>width</i>	Width of the image.
<i>height</i>	Height of the image.

Returns

int representing the average amount of red in the image.

6.4.2.2 decodeYUV420SPtoRedSum()

```
static int com.example.RUHEALTHY.demol.ImageProcessing.decodeYUV420SPtoRedSum (
    byte [] yuv420sp,
    int width,
    int height ) [static], [private]
```

Here is the caller graph for this function:

The documentation for this class was generated from the following file:

- **ImageProcessing.java**

6.5 com.example.RUHEALTHY.demo1.MainActivity ClassReference

Inherits AppCompatActivity.

Collaboration diagram for com.example.RUHEALTHY.demo1.MainActivity:

Public Member Functions

- void **LoginFunction** ()
- void **CheckEditTextStatus**()
- void **CheckFinalResult**()

Static Public Attributes

- static final String **UserEmail** = ""

Protected Member Functions

- void **onCreate** (Bundle savedInstanceState)

Package Attributes

- Button **LoginButton**
- Button **RegisterButton**
- EditText **Email**
- EditText **Password**
- String **EmailHolder**
- String **PasswordHolder**
- Boolean **EditTextEmptyHolder**
- SQLiteDatabase **sqliteDatabaseObj**
- **SQLiteHelper** **sqliteHelper**
- Cursor **cursor**
- String **TempPassword** = "NOT_FOUND"

6.5.1 Member Function Documentation

6.5.1.1 CheckEditTextStatus()

```
void com.example.RUHEALTHY.demo1.MainActivity.CheckEditTextStatus ( )
```

Here is the caller graph for this function:

6.5.1.2 CheckFinalResult()

```
void com.example.RUHEALTHY.demol.MainActivity.CheckFinalResult ( )
```

Here is the caller graph for this function:

6.5.1.3 LoginFunction()

```
void com.example.RUHEALTHY.demol.MainActivity.LoginFunction ( )
```

Here is the caller graph for this function:

6.5.1.4 onCreate()

```
void com.example.RUHEALTHY.demol.MainActivity.onCreate (
    Bundle savedInstanceState ) [protected]
```

6.5.2 Member Data Documentation

6.5.2.1 cursor

```
Cursor com.example.RUHEALTHY.demol.MainActivity.cursor [package]
```

6.5.2.2 EditTextEmptyHolder

```
Boolean com.example.RUHEALTHY.demol.MainActivity.EditTextEmptyHolder [package]
```

6.5.2.3 Email

```
EditText com.example.RUHEALTHY.demol.MainActivity.Email [package]
```

6.5.2.4 EmailHolder

```
String com.example.RUHEALTHY.demol.MainActivity.EmailHolder [package]
```

6.5.2.5 LoginButton

Button com.example.RUHEALTHY.demo1.MainActivity.LoginButton [package]

6.5.2.6 Password

EditText com.example.RUHEALTHY.demo1.MainActivity.Password [package]

6.5.2.7 PasswordHolder

String com.example.RUHEALTHY.demo1.MainActivity.PasswordHolder [package]

6.5.2.8 RegisterButton

Button com.example.RUHEALTHY.demo1.MainActivity.RegisterButton [package]

6.5.2.9 SQLiteDatabaseObj

SQLiteDatabase com.example.RUHEALTHY.demo1.MainActivity.sqliteDatabaseObj [package]

6.5.2.10 sqliteHelper

SQLiteHelper com.example.RUHEALTHY.demo1.MainActivity.sqliteHelper [package]

6.5.2.11 TempPassword

String com.example.RUHEALTHY.demo1.MainActivity.TempPassword = "NOT_FOUND" [package]

6.5.2.12 UserEmail

```
final String com.example.RUHEALTHY.demo1.MainActivity.UserEmail = "" [static]
```

The documentation for this class was generated from the following file:

- **MainActivity.java**

6.6 com.example.RUHEALTHY.demo1.ProfileActivityClass Reference

Inherits AppCompatActivity.

Collaboration diagram for com.example.RUHEALTHY.demo1.ProfileActivity:

Public Member Functions

- void **SQLiteDataBaseBuild()**
- void **InsertDataIntoSQLiteDatabase()**
- void **CheckEditTextStatus()**

Protected Member Functions

- void **onCreate** (Bundle savedInstanceState)

Package Attributes

- EditText **Age**
- EditText **Gender**
- EditText **Height**
- EditText **Weight**
- EditText **Target**
- Button **Submit**
- String **AgeHolder**
- String **GenderHolder**
- String **HeightHolder**
- String **WeightHolder**
- String **TargetHolder**
- String **EmailHolder**
- Boolean **EditTextEmptyHolder**
- SQLiteDatabase **sqliteDatabaseObj**
- String **SQLiteDataBaseQueryHolder**
- **SQLiteHelper** **sqliteHelper**

6.6.1 Member Function Documentation

6.6.1.1 CheckEditTextStatus()

```
void com.example.RUHEALTHY.demo1.ProfileActivity.CheckEditTextStatus ( )
```

Here is the caller graph for this function:

6.6.1.2 InsertDataIntoSQLiteDatabase()

```
void com.example.RUHEALTHY.demo1.ProfileActivity.InsertDataIntoSQLiteDatabase ( )
```

Here is the caller graph for this function:

6.6.1.3 onCreate()

```
void com.example.RUHEALTHY.demo1.ProfileActivity.onCreate (
    Bundle savedInstanceState ) [protected]
```

6.6.1.4 SQLiteDatabaseBuild()

```
void com.example.RUHEALTHY.demo1.ProfileActivity.SQLiteDatabaseBuild ( )
```

Here is the caller graph for this function:

6.6.2 Member Data Documentation

6.6.2.1 Age

```
EditText com.example.RUHEALTHY.demo1.ProfileActivity.Age [package]
```

6.6.2.2 AgeHolder

```
String com.example.RUHEALTHY.demo1.ProfileActivity.AgeHolder [package]
```

6.6.2.3 EditTextEmptyHolder

```
Boolean com.example.RUHEALTHY.demo1.ProfileActivity.EditTextEmptyHolder [package]
```

6.6.2.4 EmailHolder

`String com.example.RUHEALTHY.demol.ProfileActivity.EmailHolder [package]`

6.6.2.5 Gender

`EditText com.example.RUHEALTHY.demol.ProfileActivity.Gender [package]`

6.6.2.6 GenderHolder

`String com.example.RUHEALTHY.demol.ProfileActivity.GenderHolder [package]`

6.6.2.7 Height

`EditText com.example.RUHEALTHY.demol.ProfileActivity.Height [package]`

6.6.2.8 HeightHolder

`String com.example.RUHEALTHY.demol.ProfileActivity.HeightHolder [package]`

6.6.2.9 sqLiteDatabaseObj

`SQLiteDatabase com.example.RUHEALTHY.demol.ProfileActivity.sqliteDatabaseObj [package]`

6.6.2.10 SQLiteDataBaseQueryHolder

`String com.example.RUHEALTHY.demol.ProfileActivity.SQLiteDataBaseQueryHolder [package]`

6.6.2.11 sqLiteHelper

`SQLiteHelper com.example.RUHEALTHY.demol.ProfileActivity.sqliteHelper [package]`

6.6.2.12 Submit

Button com.example.RUHEALTHY.demo1.ProfileActivity.Submit [package]

6.6.2.13 Target

EditText com.example.RUHEALTHY.demo1.ProfileActivity.Target [package]

6.6.2.14 TargetHolder

String com.example.RUHEALTHY.demo1.ProfileActivity.TargetHolder [package]

6.6.2.15 Weight

EditText com.example.RUHEALTHY.demo1.ProfileActivity.Weight [package]

6.6.2.16 WeightHolder

String com.example.RUHEALTHY.demo1.ProfileActivity.WeightHolder [package]

The documentation for this class was generated from the following file:

- **ProfileActivity.java**

6.7 com.example.RUHEALTHY.demo1.RegisterActivity ClassReference

Inherits AppCompatActivity.

Collaboration diagram for com.example.RUHEALTHY.demo1.RegisterActivity:

Public Member Functions

- void **SQLiteDataBaseBuild()**
- void **SQLiteTableBuild()**
- void **InsertDataIntoSQLiteDatabase()**
- void **EmptyEditTextAfterDataInsert()**
- void **CheckEditTextStatus()**
- void **CheckingEmailAlreadyExistsOrNot()**
- void **CheckFinalResult()**

Protected Member Functions

- void **onCreate** (Bundle savedInstanceState)

Package Attributes

- EditText **Email**
- EditText **Password**
- EditText **Name**
- Button **Register**
- String **NameHolder**
- String **EmailHolder**
- String **PasswordHolder**
- Boolean **EditTextEmptyHolder**
- SQLiteDatabase **sqliteDatabaseObj**
- String **SQLiteDataBaseQueryHolder**
- **SQLiteHelper** **sqliteHelper**
- Cursor **cursor**
- String **F_Result** = "Not_Found"

6.7.1 Member Function Documentation

6.7.1.1 CheckEditTextStatus()

```
void com.example.RUHEALTHY.demol.RegisterActivity.CheckEditTextStatus ( )
```

Here is the caller graph for this function:

6.7.1.2 CheckFinalResult()

```
void com.example.RUHEALTHY.demol.RegisterActivity.CheckFinalResult ( )
```

Here is the caller graph for this function:

6.7.1.3 CheckingEmailAlreadyExistsOrNot()

```
void com.example.RUHEALTHY.demol.RegisterActivity.CheckingEmailAlreadyExistsOrNot ( )
```

Here is the caller graph for this function:

6.7.1.4 EmptyEditTextAfterDataInsert()

```
void com.example.RUHEALTHY.demol.RegisterActivity.EmptyEditTextAfterDataInsert ( )
```

Here is the caller graph for this function:

6.7.1.5 InsertDataIntoSQLiteDatabase()

```
void com.example.RUHEALTHY.demo1.RegisterActivity.InsertDataIntoSQLiteDatabase ( )
```

Here is the caller graph for this function:

6.7.1.6 onCreate()

```
void com.example.RUHEALTHY.demo1.RegisterActivity.onCreate (
    Bundle savedInstanceState ) [protected]
```

6.7.1.7 SQLiteDatabaseBuild()

```
void com.example.RUHEALTHY.demo1.RegisterActivity.SQLiteDatabaseBuild ( )
```

Here is the caller graph for this function:

6.7.1.8 SQLiteTableBuild()

```
void com.example.RUHEALTHY.demo1.RegisterActivity.SQLiteTableBuild ( )
```

Here is the caller graph for this function:

6.7.2 Member Data Documentation

6.7.2.1 cursor

```
Cursor com.example.RUHEALTHY.demo1.RegisterActivity.cursor [package]
```

6.7.2.2 EditTextEmptyHolder

```
Boolean com.example.RUHEALTHY.demo1.RegisterActivity.EditTextEmptyHolder [package]
```

6.7.2.3 Email

```
EditText com.example.RUHEALTHY.demo1.RegisterActivity.Email [package]
```

6.7.2.4 EmailHolder

```
String com.example.RUHEALTHY.demol.RegisterActivity.EmailHolder [package]
```

6.7.2.5 F_Result

```
String com.example.RUHEALTHY.demol.RegisterActivity.F_Result = "Not_Found" [package]
```

6.7.2.6 Name

```
EditText com.example.RUHEALTHY.demol.RegisterActivity.Name [package]
```

6.7.2.7 NameHolder

```
String com.example.RUHEALTHY.demol.RegisterActivity.NameHolder [package]
```

6.7.2.8 Password

```
EditText com.example.RUHEALTHY.demol.RegisterActivity.Password [package]
```

6.7.2.9 PasswordHolder

```
String com.example.RUHEALTHY.demol.RegisterActivity.PasswordHolder [package]
```

6.7.2.10 Register

```
Button com.example.RUHEALTHY.demol.RegisterActivity.Register [package]
```

6.7.2.11 sqLiteDatabaseObj

```
SQLiteDatabase com.example.RUHEALTHY.demol.RegisterActivity.sqliteDatabaseObj [package]
```

6.7.2.12 SQLiteDatabaseQueryHolder

```
String com.example.RUHEALTHY.demo1.RegisterActivity.SQLiteDatabaseQueryHolder [package]
```

6.7.2.13 sqLiteHelper

```
sqLiteHelper com.example.RUHEALTHY.demo1.RegisterActivity.sqliteHelper [package]
```

The documentation for this class was generated from the following file:

- **RegisterActivity.java**

6.8 com.example.RUHEALTHY.demo1.SensorFilter ClassReference

Static Public Member Functions

- static float **sum** (float[] array)
- static float [] **cross** (float[] arrayA, float[] arrayB)
- static float **norm** (float[] array)
- static float **dot** (float[] a, float[] b)
- static float [] **normalize** (float[] a)

Private Member Functions

- **SensorFilter ()**

6.8.1 Detailed Description

Created by aymoo on 9/19/2017.

6.8.2 Constructor & Destructor Documentation

6.8.2.1 SensorFilter()

```
com.example.RUHEALTHY.demo1.SensorFilter.SensorFilter ( ) [private]
```

6.8.3 Member Function Documentation

6.8.3.1 cross()

```
static float [ ] com.example.RUHEALTHY.demo1.SensorFilter.cross (
    float [ ] arrayA,
    float [ ] arrayB ) [static]
```

6.8.3.2 dot()

```
static float com.example.RUHEALTHY.demo1.SensorFilter.dot (
    float [ ] a,
    float [ ] b ) [static]
```

Here is the caller graph for this function:

6.8.3.3 norm()

```
static float com.example.RUHEALTHY.demo1.SensorFilter.norm (
    float [ ] array ) [static]
```

Here is the caller graph for this function:

6.8.3.4 normalize()

```
static float [ ] com.example.RUHEALTHY.demo1.SensorFilter.normalize (
    float [ ] a ) [static]
```

6.8.3.5 sum()

```
static float com.example.RUHEALTHY.demo1.SensorFilter.sum (
    float [ ] array ) [static]
```

Here is the caller graph for this function:

The documentation for this class was generated from the following file:

- **SensorFilter.java**

6.9 com.example.RUHEALTHY.demo1.SQLiteHelper ClassReference

Inherits SQLiteOpenHelper.

Public Member Functions

- **SQLiteHelper** (Context context)
- void **onCreate** (SQLiteDatabase database)
- void **onUpgrade** (SQLiteDatabase db, int oldVersion, int newVersion)

Static Public Attributes

- static final String **TABLE_NAME**="usertable"
- static final String **Table_Column_1_Name**="name"
- static final String **Table_Column_2_Email**="email"
- static final String **Table_Column_3_Password**="password"
- static final String **Table_Column_4_Age**="age"
- static final String **Table_Column_5_Gender**="gender"
- static final String **Table_Column_6_Height**="height"
- static final String **Table_Column_7_Weight**="weight"
- static final String **Table_Column_8_Target**="target"

Static Package Attributes

- static String **DATABASE_NAME**="UserDataBase"

6.9.1 Detailed Description

Created by RUHEALTHY on 10/27/2017.

6.9.2 Constructor & Destructor Documentation

6.9.2.1 SQLiteHelper()

```
com.example.RUHEALTHY.demo1.SQLiteHelper.SQLiteHelper (
    Context context )
```

6.9.3 Member Function Documentation

6.9.3.1 onCreate()

```
void com.example.RUHEALTHY.demo1.SQLiteHelper.onCreate (
    SQLiteDatabase database )
```

Here is the caller graph for this function:

6.9.3.2 onUpgrade()

```
void      com.example.RUHEALTHY.demo1.SQLiteHelper.onUpgrade      (  
        SQLiteDatabase db,  
        int oldVersion,  
        int newVersion )
```

6.9.4 Member Data Documentation

6.9.4.1 DATABASE_NAME

```
String com.example.RUHEALTHY.demo1.SQLiteHelper.DATABASE_NAME = "UserDataBase"  [static], [package]
```

6.9.4.2 Table_Column_1_Name

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_1_Name = "name" [static]
```

6.9.4.3 Table_Column_2_Email

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_2_Email = "email" [static]
```

6.9.4.4 Table_Column_3_Password

```
final      String      com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_3_Password      = "password"  
[static]
```

6.9.4.5 Table_Column_4_Age

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_4_Age = "age" [static]
```

6.9.4.6 Table_Column_5_Gender

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_5_Gender = "gender" [static]
```

6.9.4.7 Table_Column_6_Height

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_6_Height ="height" [static]
```

6.9.4.8 Table_Column_7_Weight

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_7_Weight ="weight" [static]
```

6.9.4.9 Table_Column_8_Target

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.Table_Column_8_Target ="target" [static]
```

6.9.4.10 TABLE_NAME

```
final String com.example.RUHEALTHY.demo1.SQLiteHelper.TABLE_NAME ="usertable" [static]
```

The documentation for this class was generated from the following file:

- **SQLiteHelper.java**

6.10 com.example.RUHEALTHY.demo1.StepCounter ClassReference

Inheritance diagram for com.example.RUHEALTHY.demo1.StepCounter:

Collaboration diagram for

com.example.RUHEALTHY.demo1.StepCounter:

Public Member Functions

- void **onSensorChanged** (SensorEventevent)
- void **onAccuracyChanged** (Sensor sensor, int accuracy)
- void **step** (longtimeNs)

Protected Member Functions

- void **onCreate** (BundlesavedInstanceState)

Private Attributes

- TextView **TvMiles**
- TextView **TvSteps**
- TextView **TvCalories**
- **StepDetector** **simpleStepDetector**
- SensorManager **sensorManager**
- Sensor **accel**
- double **numSteps**
- double **miles**
- double **calories**

Static Private Attributes

- static final String **TEXT_NUM_STEPS** = "STEP"
- static final String **TEXT_Miles_STEPS** = "MILE"
- static final String **TEXT_Miles_Calorie** = "CALORIES PER MILE "

6.10.1 Member Function Documentation

6.10.1.1 onAccuracyChanged()

```
void    com.example.RUHEALTHY.demo1.StepCounter.onAccuracyChanged    (
        Sensor sensor,
        int accuracy )
```

6.10.1.2 onCreate()

```
void    com.example.RUHEALTHY.demo1.StepCounter.onCreate    (
        Bundle savedInstanceState ) [protected]
```

6.10.1.3 onSensorChanged()

```
void    com.example.RUHEALTHY.demo1.StepCounter.onSensorChanged    (
        SensorEvent event )
```

6.10.1.4 step()

```
void    com.example.RUHEALTHY.demo1.StepCounter.step    (
        long timeNs )
```

Implements **com.example.RUHEALTHY.demo1.StepListener** (p. ??).

6.10.2 Member Data Documentation

6.10.2.1 accel

Sensor com.example.RUHEALTHY.demo1.StepCounter.accel [private]

6.10.2.2 calories

double com.example.RUHEALTHY.demo1.StepCounter.calories [private]

6.10.2.3 miles

double com.example.RUHEALTHY.demo1.StepCounter.miles [private]

6.10.2.4 numSteps

double com.example.RUHEALTHY.demo1.StepCounter.numSteps [private]

6.10.2.5 sensorManager

SensorManager com.example.RUHEALTHY.demo1.StepCounter.sensorManager [private]

6.10.2.6 simpleStepDetector

StepDetector com.example.RUHEALTHY.demo1.StepCounter.simpleStepDetector [private]

6.10.2.7 TEXT_Miles_Calorie

final String com.example.RUHEALTHY.demo1.StepCounter.TEXT_Miles_Calorie = "CALORIES PER MILE "
[static], [private]

6.10.2.8 TEXT_Miles_STEPS

```
final String com.example.RUHEALTHY.demo1.StepCounter.TEXT_Miles_STEPS = "MILE " [static],  
[private]
```

6.10.2.9 TEXT_NUM_STEPS

```
final String com.example.RUHEALTHY.demo1.StepCounter.TEXT_NUM_STEPS = "STEP " [static],[private]
```

6.10.2.10 TvCalories

```
TextView com.example.RUHEALTHY.demo1.StepCounter.TvCalories [private]
```

6.10.2.11 TvMiles

```
TextView com.example.RUHEALTHY.demo1.StepCounter.TvMiles [private]
```

6.10.2.12 TvSteps

```
TextView com.example.RUHEALTHY.demo1.StepCounter.TvSteps [private]
```

The documentation for this class was generated from the following file:

- **StepCounter.java**

6.11 com.example.RUHEALTHY.demo1.StepDetector ClassReference

Collaboration diagram for com.example.RUHEALTHY.demo1.StepDetector:

Public Member Functions

- void **registerListener** (**StepListener** listener)
- void **updateAccel** (long timeNs, float x, float y, float z)

Private Attributes

- int **accelRingCounter** = 0
- float [] **accelRingX** = new float[**ACCEL_RING_SIZE**]
- float [] **accelRingY** = new float[**ACCEL_RING_SIZE**]
- float [] **accelRingZ** = new float[**ACCEL_RING_SIZE**]
- int **velRingCounter** = 0
- float [] **velRing** = new float[**VEL_RING_SIZE**]
- long **lastStepTimeNs** = 0
- float **oldVelocityEstimate** = 0
- **StepListener** listener

Static Private Attributes

- static final int **ACCEL_RING_SIZE** = 50
- static final int **VEL_RING_SIZE** = 10
- static final float **STEP_THRESHOLD** = 50f
- static final int **STEP_DELAY_NS** = 250000000

6.11.1 Detailed Description

Created by aymoo on 9/19/2017.

6.11.2 Member Function Documentation

6.11.2.1 registerListener()

```
void com.example.RUHEALTHY.demo1.StepDetector.registerListener (
    StepListener listener )
```

Here is the caller graph for this function:

6.11.2.2 updateAccel()

```
void com.example.RUHEALTHY.demo1.StepDetector.updateAccel (
    long timeNs,
    float x,
    float y,
    float z )
```

Here is the caller graph for this function:

6.11.3 Member Data Documentation

6.11.3.1 ACCEL_RING_SIZE

```
final int com.example.RUHEALTHY.demol.StepDetector.ACCEL_RING_SIZE = 50 [static],[private]
```

6.11.3.2 accelRingCounter

```
int com.example.RUHEALTHY.demol.StepDetector.accelRingCounter =0 [private]
```

6.11.3.3 accelRingX

```
float [ ] com.example.RUHEALTHY.demol.StepDetector.accelRingX = new float[ ACCEL_RING_SIZE]  
[private]
```

6.11.3.4 accelRingY

```
float [ ] com.example.RUHEALTHY.demol.StepDetector.accelRingY = new float[ ACCEL_RING_SIZE]  
[private]
```

6.11.3.5 accelRingZ

```
float [ ] com.example.RUHEALTHY.demol.StepDetector.accelRingZ = new float[ ACCEL_RING_SIZE]  
[private]
```

6.11.3.6 lastStepTimeNs

```
long com.example.RUHEALTHY.demol.StepDetector.lastStepTimeNs = 0 [private]
```

6.11.3.7 listener

```
StepListener com.example.RUHEALTHY.demol.StepDetector.listener [private]
```

6.11.3.8 oldVelocityEstimate

```
float com.example.RUHEALTHY.demo1.StepDetector.oldVelocityEstimate = 0 [private]
```

6.11.3.9 STEP_DELAY_NS

```
final int com.example.RUHEALTHY.demo1.StepDetector.STEP_DELAY_NS = 250000000 [static],[private]
```

6.11.3.10 STEP_THRESHOLD

```
final float com.example.RUHEALTHY.demo1.StepDetector.STEP_THRESHOLD = 50f [static],[private]
```

6.11.3.11 VEL_RING_SIZE

```
final int com.example.RUHEALTHY.demo1.StepDetector.VEL_RING_SIZE = 10 [static],[private]
```

6.11.3.12 velRing

```
float [ ] com.example.RUHEALTHY.demo1.StepDetector.velRing = new float[ VEL_RING_SIZE] [private]
```

6.11.3.13 velRingCounter

```
int com.example.RUHEALTHY.demo1.StepDetector.velRingCounter = 0 [private]
```

The documentation for this class was generated from the following file:

- **StepDetector.java**

6.12 com.example.RUHEALTHY.demo1.StepListener InterfaceReference

Inheritance diagram for com.example.RUHEALTHY.demo1.StepListener:

Public Member Functions

- void **step** (long timeNs)

6.12.1 Detailed Description

Created by aymoo on 9/19/2017.

6.12.2 Member Function Documentation

6.12.2.1 step()

```
void com.example.RUHEALTHY.demo1.StepListener.step (
    long timeNs )
```

Implemented in **com.example.RUHEALTHY.demo1.StepCounter** (p.

??). Here is the caller graph for this function:

The documentation for this interface was generated from the following file:

- **StepListener.java**

6.13 com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPEEnumReference

Public Attributes

- **GREEN**
- **RED**

6.13.1 Member Data Documentation

6.13.1.1 GREEN

```
com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE.GREEN
```

6.13.1.2 RED

```
com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE.RED
```

The documentation for this enum was generated from the following file:

- **HeartRateMonitor.java**

Chapter 7

File Documentation

7.1 DashboardActivity.java File Reference

Classes

- class **com.example.RUHEALTHY.demo1.DashboardActivity**

Packages

- package **com.example.RUHEALTHY.demo1**

7.2 HeartbeatView.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.HeartbeatView**

Packages

- package **com.example.RUHEALTHY.demo1**

7.3 HeartRateMonitor.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.HeartRateMonitor**
- enum **com.example.RUHEALTHY.demo1.HeartRateMonitor.TYPE**

Packages

- package **com.example.RUHEALTHY.demo1**

7.4 ImageProcessing.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.ImageProcessing**

Packages

- package **com.example.RUHEALTHY.demo1**

7.5 MainActivity.java File Reference

Classes

- class **com.example.RUHEALTHY.demo1.MainActivity**

Packages

- package **com.example.RUHEALTHY.demo1**

7.6 ProfileActivity.java File Reference

Classes

- class **com.example.RUHEALTHY.demo1.ProfileActivity**

Packages

- package **com.example.RUHEALTHY.demo1**

7.7 RegisterActivity.java File Reference

Classes

- class **com.example.RUHEALTHY.demo1.RegisterActivity**

Packages

- package **com.example.RUHEALTHY.demo1**

7.8 SensorFilter.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.SensorFilter**

Packages

- package **com.example.RUHEALTHY.demo1**

7.9 SQLiteHelper.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.SQLiteHelper**

Packages

- package **com.example.RUHEALTHY.demo1**

7.10 StepCounter.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.StepCounter**

Packages

- package **com.example.RUHEALTHY.demo1**

7.11 StepDetector.java FileReference

Classes

- class **com.example.RUHEALTHY.demo1.StepDetector**

Packages

- package **com.example.RUHEALTHY.demo1**

7.12 StepListener.java FileReference

Classes

- interface **com.example.RUHEALTHY.demo1.StepListener**

Packages

- package **com.example.RUHEALTHY.demo1**

Index

ACCEL_RING_SIZE	com::example::RUHEALTHY::demo1::StepDetector, 39
accel	com::example::RUHEALTHY::demo1::StepCounter, 37
accelRingCounter	com::example::RUHEALTHY::demo1::StepDetector, 40
accelRingX	com::example::RUHEALTHY::demo1::StepDetector, 40
accelRingY	com::example::RUHEALTHY::demo1::StepDetector, 40
accelRingZ	com::example::RUHEALTHY::demo1::StepDetector, 40
Age	com::example::RUHEALTHY::demo1::ProfileActivity, 25
AgeHolder	com::example::RUHEALTHY::demo1::ProfileActivity, 25
averageArray	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
averageArraySize	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
averageIndex	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
beats	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
beatsArray	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
beatsArraySize	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
beatsIndex	com::example::RUHEALTHY::demo1::HeartRateMonitor, 17
calories	com::example::RUHEALTHY::demo1::StepCounter, 37

beatsArraySize, 17

camera

com::example::RUHEALTHY::demo1::HeartRate
 ← 'Monitor, 18

CheckEditTextStatus

com::example::RUHEALTHY::demo1::MainActivi
 ty, 21
 com::example::RUHEALTHY::demo1::ProfileActi
 vity,
 24
 com::example::RUHEALTHY::demo1::Register← '
 Activity, 28

CheckFinalResult

com::example::RUHEALTHY::demo1::MainActivi
 ty, 21
 com::example::RUHEALTHY::demo1::Register← '
 Activity, 28

CheckingEmailAlreadyExistsOrNot

com::example::RUHEALTHY::demo1::Register← '
 Activity, 28

com, 9

com.example, 9

com.example.RUHEALTHY, 9

com.example.RUHEALTHY.demo1, 9

com.example.RUHEALTHY.demo1.DashboardActivit
 y, 11

com.example.RUHEALTHY.demo1.HeartRateMonitor
 , 15

com.example.RUHEALTHY.demo1.HeartRateMonitor
 .T← '

YPE,
 42

com.example.RUHEALTHY.demo1.HeartbeatView,
 12

com.example.RUHEALTHY.demo1.ImageProcessing,
 19

com.example.RUHEALTHY.demo1.MainActivity, 21

com.example.RUHEALTHY.demo1.ProfileActivity, 24

com.example.RUHEALTHY.demo1.RegisterActivity,
 27

com.example.RUHEALTHY.demo1.SQLiteHelper, 32

com.example.RUHEALTHY.demo1.SensorFilter, 31

com.example.RUHEALTHY.demo1.StepCounter, 35

com.example.RUHEALTHY.demo1.StepDetector, 38

com.example.RUHEALTHY.demo1.StepListener, 41

com::example::RUHEALTHY::demo1::DashboardActi
 vity

Email, 11

EmailHolder, 12

HRM, 12

LogOUT, 12

onCreate, 11

Profile, 12

Step, 12

com::example::RUHEALTHY::demo1::HeartRateMonitor

averageArray, 17

averageArraySize, 17

averageIndex, 17

beats, 17

beatsArray, 17

- beatsIndex, 17
- camera, 18
- currentType, 18
- getCurrent, 16
- getSmallestPreviewSize, 16
- image, 18
- onConfigurationChanged, 16
- onCreate, 16
- onPause, 16
- onResume, 16
- preview, 18
- previewCallback, 18
- previewHolder, 18
- processing, 18
- startTime, 19
- surfaceCallback, 19
- TAG, 19
- text, 19
- wakeLock, 19
- com::example::RUHEALTHY::demo1::HeartRateMonitor
 - TYPE
 - GREEN, 42
 - RED, 42
- com::example::RUHEALTHY::demo1::HeartbeatView
 - greenBitmap, 14
 - HeartbeatView, 13
 - matrix, 14
 - onDraw, 13
 - onMeasure, 14
 - paint, 14
 - parentHeight, 14
 - parentWidth, 14
 - redBitmap, 14
- com::example::RUHEALTHY::demo1::ImageProcessing
 - decodeYUV420SPtoRedAvg, 20
 - decodeYUV420SPtoRedSum, 20
- com::example::RUHEALTHY::demo1::MainActivity
 - CheckEditTextStatus, 21
 - CheckFinalResult, 21
 - cursor, 22
 - EditTextEmptyHolder, 22
 - Email, 22
 - EmailHolder, 22
 - LogInButton, 22
 - LoginFunction, 22
 - onCreate, 22
 - Password, 23
 - PasswordHolder, 23
 - RegisterButton, 23
 - sqliteDatabaseObj, 23
 - sqliteHelper, 23
 - TempPassword, 23
 - UserEmail, 23
- com::example::RUHEALTHY::demo1::ProfileActivity
 - Age, 25
 - AgeHolder, 25
 - CheckEditTextStatus, 24
- EditTextEmptyHolder, 25

- EmailHolder, 25
- Gender, 26
- GenderHolder, 26
- Height, 26
- HeightHolder, 26
- InsertDataIntoSQLiteDatabase, 25
- onCreate, 25
- SQLiteDataBaseBuild, 25
- SQLiteDataBaseQueryHolder, 26
- sqliteDatabaseObj, 26
- sqliteHelper, 26
- Submit, 26
- Target, 27
- TargetHolder, 27
- Weight, 27
- WeightHolder, 27
- com::example::RUHEALTHY::demo1::RegisterActivity
 - CheckEditTextStatus, 28
 - CheckFinalResult, 28
 - CheckingEmailAlreadyExistsOrNot, 28
 - cursor, 29
 - EditTextEmptyHolder, 29
 - Email, 29
 - EmailHolder, 29
 - EmptyEditTextAfterDataInsert, 28
 - F_Result, 30
 - InsertDataIntoSQLiteDatabase, 28
 - Name, 30
 - NameHolder, 30
 - onCreate, 29
 - Password, 30
 - PasswordHolder, 30
 - Register, 30
 - SQLiteDataBaseBuild, 29
 - SQLiteDataBaseQueryHolder, 30
 - SQLiteTableBuild, 29
 - sqliteDatabaseObj, 30
 - sqliteHelper, 31
- com::example::RUHEALTHY::demo1::SQLiteHelper
 - DATABASE_NAME, 34
 - onCreate, 33
 - onUpgrade, 33
 - SQLiteHelper, 33
 - TABLE_NAME, 35
 - Table_Column_1_Name, 34
 - Table_Column_2_Email, 34
 - Table_Column_3_Password, 34
 - Table_Column_4_Age, 34
 - Table_Column_5_Gender, 34
 - Table_Column_6_Height, 34
 - Table_Column_7_Weight, 35
 - Table_Column_8_Target, 35
- com::example::RUHEALTHY::demo1::SensorFilter
 - cross, 31
 - dot, 32
 - norm, 32
 - normalize, 32
 - SensorFilter, 31

sum, 32 dot
 com::example::RUHEALTHY::demo1::StepCounter
 accel, 37
 calories, 37
 miles, 37
 numSteps, 37
 onAccuracyChanged, 36
 onCreate, 36
 onSensorChanged, 36
 sensorManager, 37
 simpleStepDetector, 37
 step, 36
 TEXT_Miles_Calorie, 37
 TEXT_Miles_STEPS, 37
 TEXT_NUM_STEPS, 38
 TvCalories, 38
 TvMiles, 38
 TvSteps, 38
 com::example::RUHEALTHY::demo1::StepDetector
 ACCEL_RING_SIZE, 39
 accelRingCounter, 40
 accelRingX, 40
 accelRingY, 40
 accelRingZ, 40
 lastStepTimeNs, 40
 listener, 40
 oldVelocityEstimate, 40
 registerListener, 39
 STEP_DELAY_NS, 41
 STEP_THRESHOLD, 41
 updateAccel, 39
 VEL_RING_SIZE, 41
 velRing, 41
 velRingCounter, 41
 com::example::RUHEALTHY::demo1::StepListener
 step, 42
 cross
 com::example::RUHEALTHY::demo1::SensorFilter, 31
 currentType
 com::example::RUHEALTHY::demo1::HeartRateMonitor, 18
 Monitor, 18
 cursor
 com::example::RUHEALTHY::demo1::MainActivity, 22
 com::example::RUHEALTHY::demo1::RegisterActivity, 29

 DATABASE_NAME
 com::example::RUHEALTHY::demo1::SQLiteHelper, 34
 DashboardActivity.java, 43
 decodeYUV420SPtoRedAvg
 com::example::RUHEALTHY::demo1::ImageProcessing, 20
 decodeYUV420SPtoRedSum
 com::example::RUHEALTHY::demo1::ImageProcessing, 20

- com::example::RUHEALTHY::demo1::SensorFilter, 32
- EditTextEmptyHolder
 - com::example::RUHEALTHY::demo1::MainActivity, 22
 - com::example::RUHEALTHY::demo1::ProfileActivity, 25
 - com::example::RUHEALTHY::demo1::Register←'Activity, 29
- Email
 - com::example::RUHEALTHY::demo1::Dashboard←'Activity, 11
 - com::example::RUHEALTHY::demo1::MainActivity, 22
 - com::example::RUHEALTHY::demo1::Register←'Activity, 29
- 29 EmailHolder
 - com::example::RUHEALTHY::demo1::Dashboard←'Activity, 12
 - com::example::RUHEALTHY::demo1::MainActivity, 22
 - com::example::RUHEALTHY::demo1::ProfileActivity, 25
 - com::example::RUHEALTHY::demo1::Register←'Activity, 29
- EmptyEditTextAfterDataInsert
 - com::example::RUHEALTHY::demo1::Register←'Activity, 28
- F_Result
 - com::example::RUHEALTHY::demo1::Register←'Activity, 30
- GREEN
 - com::example::RUHEALTHY::demo1::HeartRate←'Monitor::TYPE, 42
- Gender
 - com::example::RUHEALTHY::demo1::ProfileActivity, 26
- GenderHolder
 - com::example::RUHEALTHY::demo1::ProfileActivity, 26
- getCurrent
 - com::example::RUHEALTHY::demo1::HeartRate←'Monitor, 16
- getSmallestPreviewSize
 - com::example::RUHEALTHY::demo1::HeartRate←'Monitor, 16
- 16 greenBitmap
 - com::example::RUHEALTHY::demo1::HeartbeatView, 14
- HRM
 - com::example::RUHEALTHY::demo1::Dashboard←'Activity, 12
 - HeartRateMonitor.java, 43
 - HeartbeatView
 - com::example::RUHEALTHY::demo1::HeartbeatView, 13
 - HeartbeatView.java, 43
 - Height

- com::example::RUHEALTHY::demo1::ProfileActivi
ty, 26
- HeightHolder
 - com::example::RUHEALTHY::demo1::ProfileActivi
ty, 26
- image
 - com::example::RUHEALTHY::demo1::HeartRate←
' Monitor, 18
- ImageProcessing.java, 44
- InsertDataIntoSQLiteDatabase
 - com::example::RUHEALTHY::demo1::ProfileActivi
ty, 25
 - com::example::RUHEALTHY::demo1::Register←
Activity, 28
- lastStepTimeNs
 - com::example::RUHEALTHY::demo1::StepDetect
or, 40
- listener
 - com::example::RUHEALTHY::demo1::StepDetect
or, 40
- LogInButton
 - com::example::RUHEALTHY::demo1::MainActivity
, 22
- LogOUT
 - com::example::RUHEALTHY::demo1::Dashboard
←' Activity, 12
- LoginFunction
 - com::example::RUHEALTHY::demo1::MainActivity
, 22
- MainActivity.java, 44
- matrix
 - com::example::RUHEALTHY::demo1::HeartbeatVie
w, 14
- miles
 - com::example::RUHEALTHY::demo1::StepCounte
r, 37
- Name
 - com::example::RUHEALTHY::demo1::Register←
Activity, 30
- NameHolder
 - com::example::RUHEALTHY::demo1::Register←
Activity, 30
- norm
 - com::example::RUHEALTHY::demo1::SensorFilde
r, 32
- normalize
 - com::example::RUHEALTHY::demo1::SensorFilde
r, 32
- numSteps
 - com::example::RUHEALTHY::demo1::StepCounte
r, 37
- oldVelocityEstimate
 - com::example::RUHEALTHY::demo1::StepDetec
tor,
- onConfigurationChanged
 - com::example::RUHEALTHY::demo1::HeartRate←
Monitor, 16
- onCreate
 - com::example::RUHEALTHY::demo1::Dashboard←
Activity, 11
 - com::example::RUHEALTHY::demo1::HeartRate←
Monitor, 16
 - com::example::RUHEALTHY::demo1::MainActivity,
22
 - com::example::RUHEALTHY::demo1::ProfileActivity
, 25
 - com::example::RUHEALTHY::demo1::Register←
Activity, 29
 - com::example::RUHEALTHY::demo1::SQLiteHelper
, 33
 - com::example::RUHEALTHY::demo1::StepCounte
r, 36
- onDraw
 - com::example::RUHEALTHY::demo1::HeartbeatVi
ew, 13
- onMeasure
 - com::example::RUHEALTHY::demo1::HeartbeatVi
ew, 14
- onPause
 - com::example::RUHEALTHY::demo1::HeartRate←
Monitor, 16
- onResume
 - com::example::RUHEALTHY::demo1::HeartRate←
Monitor, 16
- onSensorChanged
 - com::example::RUHEALTHY::demo1::StepCounte
r, 36
- onUpgrade
 - com::example::RUHEALTHY::demo1::SQLiteHelp
er, 33
- paint
 - com::example::RUHEALTHY::demo1::HeartbeatVi
ew, 14
- parentHeight
 - com::example::RUHEALTHY::demo1::HeartbeatVi
ew, 14
- parentWidth
 - com::example::RUHEALTHY::demo1::HeartbeatVi
ew, 14
- Password
 - com::example::RUHEALTHY::demo1::MainActivity,
23
 - com::example::RUHEALTHY::demo1::Register←
Activity, 30
- PasswordHolder
 - com::example::RUHEALTHY::demo1::MainActivity,
23
 - com::example::RUHEALTHY::demo1::Register←
Activity, 30
- onAccuracyChanged
 - com::example::RUHEALTHY::demo1::StepCounter,
36


```
preview
  com::example::RUHEALTHY::demo1::HeartRate
  ←'
    Monitor
,      18
previewCallback
  com::example::RUHEALTHY::demo1::HeartRate
  ←' Monitor, 18
```

previewHolder
 com::example::RUHEALTHY::demo1::HeartRate←
 Monitor, 18
 processing
 com::example::RUHEALTHY::demo1::HeartRate←
 Monitor, 18
 Profile
 com::example::RUHEALTHY::demo1::Dashboard
 ← Activity, 12
 ProfileActivity.java, 44

 RED
 com::example::RUHEALTHY::demo1::HeartRate←
 Monitor::TYPE, 42
 redBitmap
 com::example::RUHEALTHY::demo1::HeartbeatVie
 w, 14
 Register
 com::example::RUHEALTHY::demo1::Register←
 Activity, 30
 RegisterActivity.java, 44
 RegisterButton
 com::example::RUHEALTHY::demo1::MainActivity
 , 23 registerListener
 com::example::RUHEALTHY::demo1::StepDetect
 or, 39

 SQLiteDatabaseBuild
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty,
 25
 com::example::RUHEALTHY::demo1::Register←
 Activity, 29
 SQLiteDatabaseQueryBuilder
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty,
 26
 com::example::RUHEALTHY::demo1::Register←
 Activity, 30
 SQLiteHelper
 com::example::RUHEALTHY::demo1::SQLiteHelp
 er, 33
 SQLiteHelper.java, 45
 SQLiteTableBuild
 com::example::RUHEALTHY::demo1::Register←
 Activity, 29
 STEP_DELAY_NS
 com::example::RUHEALTHY::demo1::StepDetect
 or, 41
 STEP_THRESHOLD
 com::example::RUHEALTHY::demo1::StepDetect
 or, 41
 SensorFilter
 com::example::RUHEALTHY::demo1::SensorFilde
 r, 31
 SensorFilter.java, 45
 sensorManager
 com::example::RUHEALTHY::demo1::StepCounter,
 37
 simpleStepDetector

com::example::RUHEALTHY::demo1::StepCounter, 37
 SQLiteDatabaseObj
 com::example::RUHEALTHY::demo1::MainActivity, 23
 com::example::RUHEALTHY::demo1::ProfileActivity, 26
 com::example::RUHEALTHY::demo1::Register←'Activity',
 30 SQLiteHelper
 com::example::RUHEALTHY::demo1::MainActivity, 23
 com::example::RUHEALTHY::demo1::ProfileActivity, 26
 com::example::RUHEALTHY::demo1::Register←'Activity', 31
 startTime
 com::example::RUHEALTHY::demo1::HeartRate←'Monitor', 19
 Step
 com::example::RUHEALTHY::demo1::Dashboard←'Activity', 12
 step
 com::example::RUHEALTHY::demo1::StepCounter, 36
 com::example::RUHEALTHY::demo1::StepListener, 42
 StepCounter.java, 45
 StepDetector.java, 45
 StepListener.java, 45
 Submit
 com::example::RUHEALTHY::demo1::ProfileActivity, 26
 sum
 com::example::RUHEALTHY::demo1::SensorFilter, 32
 surfaceCallback
 com::example::RUHEALTHY::demo1::HeartRate←'Monitor', 19
 TABLE_NAME
 com::example::RUHEALTHY::demo1::SQLiteHelper, 35
 TAG
 com::example::RUHEALTHY::demo1::HeartRate←'Monitor', 19
 TEXT_Miles_Calorie
 com::example::RUHEALTHY::demo1::StepCounter, 37
 TEXT_Miles_STEPS
 com::example::RUHEALTHY::demo1::StepCounter, 37
 TEXT_NUM_STEPS
 com::example::RUHEALTHY::demo1::StepCounter, 38

Table_Column_1_Name
 com::example::RUHEALTHY::demo1::SQLiteHelper, 34
 Table_Column_2_Email
 com::example::RUHEALTHY::demo1::SQLiteHelper, 34

Table_Column_3_Password
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 34
 Table_Column_4_Age
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 34
 Table_Column_5_Gender
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 34
 Table_Column_6_Height
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 34
 Table_Column_7_Weight
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 35
 Table_Column_8_Target
 com::example::RUHEALTHY::demo1::SQLiteHel
 per, 35
 Target
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty, 27
 TargetHolder
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty, 27
 TempPassword
 com::example::RUHEALTHY::demo1::MainActivity, 23
 text
 com::example::RUHEALTHY::demo1::HeartRateMoni
 tor, 19
 TvCalories
 com::example::RUHEALTHY::demo1::StepCounte
 r, 38
 TvMiles
 com::example::RUHEALTHY::demo1::StepCounte
 r, 38
 TvSteps
 com::example::RUHEALTHY::demo1::StepCounte
 r, 38
 updateAccel
 com::example::RUHEALTHY::demo1::StepDetect
 or, 39
 UserEmail
 com::example::RUHEALTHY::demo1::MainActivity, 23
 VEL_RING_SIZE
 com::example::RUHEALTHY::demo1::StepDetect
 or, 41
 velRing
 com::example::RUHEALTHY::demo1::StepDetect
 or, 41
 velRingCounter
 com::example::RUHEALTHY::demo1::StepDetect
 or, 41
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty, 27
 WeightHolder
 com::example::RUHEALTHY::demo1::ProfileActivi
 ty, 27